

“1.00 g of a Compound” Lab Briefing

- A **double replacement** reaction of **two (aq)** solutions will often produce a **precipitate (s)**
 - **A:** The **solid** precipitate can then be easily isolated through a **filtration process**.
 - **B:** Some ionic compounds contain **water molecules (hydrates)** that are trapped within their crystal structure and therefore must be **dried out**.
 - **LAB OBJECTIVE:** To produce **EXACTLY 1.00 grams** of desired product (**precipitate**) using solubility rules of a double replacement reaction:
 - **Pre-Lab Step 1:** Write a **balanced chemical equation** for your given reaction.
 - **Pre-Lab Step 2:** Assign states of matter (**aq or s**) to **ALL** reactants and products in the chemical equation using **solubility rules**.
 - **Pre-Lab Step 3:** Calculate the **THEORETICAL MASS (g)** of **BOTH REACTANTS** needed to exactly produce the desired 1.00 grams of the **solid precipitate**.
 - **A:** 1.00 g of solid precipitate is the **GIVEN** → Use **MASS → MASS** stoichiometry to determine required **mass** of **BOTH REACTANTS**. Therefore, **two (2) mass → mass** stoichiometry calculation set-up required.
 - **B:** If **ANY** reactant is a **HYDRATE**, you **MUST** take the **hydrate mass** into account when determining its **total MOLAR MASS**.
 - **C:** Number of **water molecules MUST** also be written on the product side to correctly balance the equation.
 - Ex: $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}_{(aq)} + \text{Na}_2\text{CO}_{3(aq)} \rightarrow \text{MgCO}_3(s) + \text{Na}_2\text{SO}_{4(aq)}$
 - $1 \text{MgSO}_4 \cdot 7\text{H}_2\text{O}_{(aq)} + 1 \text{Na}_2\text{CO}_{3(aq)} \rightarrow 1 \text{MgCO}_3(s) + 1 \text{Na}_2\text{SO}_{4(aq)} + 7\text{H}_2\text{O}$
 - $\text{MgSO}_4 \cdot 7\text{H}_2\text{O} \rightarrow 120.37 \text{ g/mol } \text{MgSO}_4 + 126.112 \text{ g/mol } 7\text{H}_2\text{O} = 246.4822 \text{ g/mol}$
 - **Pre-Lab Step 4:** Watch the provided **VIDEO** for performing the proper filtration method/process.
-

- **Actual Lab Step 1:** Once the **THEORETICAL MASS** of both reactants are **CALCULATED**, dissolve both reactants in **SEPARATE** beakers using **25mL** of **DISTILLED** water.
 - **Actual Lab Step 2:** The two (2) **DISSOLVED** reactant solutions are then **MIXED** together in ***ONE of the two*** beakers (or into a third beaker) to produce the **solid precipitate**.
 - **Actual Lab Step 3:** **RECOVER** precipitate through **FILTRATION** method using a piece of **PRE-MASSED filter paper**.
 - **Actual Lab Step 4:** **CAREFULLY** remove **ALL** precipitate **AND** filter paper from funnel and place on **DRY** towel.
 - **Actual Lab Step 5:** **MASS** dry precipitate **AND** filter paper together. **SUBTRACT** the total mass (precipitate + filter paper) from the original mass of filter paper to get just the mass of the precipitate.
-

- **AFTER Lab Step 1:** Calculate **% YIELD** to determine how much of precipitate was **ACTUALLY** produced or recovered.
 - How close is *your* **ACTUAL** mass to the **1.00 g** **THEORETICAL** mass of precipitate?

- **% Yield** =
$$\frac{\text{Actual Yield of Precipitate}}{\text{Theoretical Yield of Precipitate}} \times 100$$

- **AFTER Lab Step 2:** Calculate **% ERROR** to determine how much of precipitate was **NOT** produced or recovered.

- **% Error** =
$$\frac{[\text{Actual} - \text{Theoretical}]}{\text{Theoretical}} \times 100$$